

# VivoForce Force sensitive hook for epiretinal membrane peeling in eye surgery

Project management: **Dr Charles Baur** [charles.baur@epfl.ch](mailto:charles.baur@epfl.ch)

Start date 1.09.2014

EPFL – STI – IMT – INSTANT-LAB


Duration 2 years

## Project goal


This project addresses the **design, construction and evaluation of a peeling hook with force measurement capability for in-vivo intra-ocular vitreoretinal surgery**. The force sensor consists of a miniature multi-degree-of-freedom flexure where deformations induced by contact forces are measured using optical fiber white light interferometry. This instrument will be used for epiretinal membrane peeling procedures and should then lead to the creation of a new generation of force sensitive surgical tools.

## Medical need

Forces exerted onto retinal structures may generate irreversible visual impairment if they exceed specific thresholds. Since these forces are too small to be detected by the tactile senses of the surgeon's hand, they are nowadays estimated based solely on the visual feedback through the microscope. The novel instrument will thus improve patient safety by providing an objective measurement of the exerted forces.



## Proposed solution





### In-vivo force sensing surgical instrument

- Force range: 0-15mN
- Measuring resolution: 0.1mN
- Tool diameter: 0.6 – 0.9 mm
- Biocompatible, adapted for liquid environment
- Immune to electric and magnetic noise

### Portable interface

- Real-time force value monitoring and chart display
- Audio feedback indicating safety threshold

## Force measurement using white light interferometry



## Forseen benefits

### Safer vitreoretinal surgery

- Real-time feedback of peeling forces
- Information on forces below human perception
- Quantification of the tissue properties

### Improved training of surgeons

- Correlation of surgical outcome with exerted forces
- Faster learning curve
- Standardization: correlate forces with surgical parameters

### Overall clinical impact

- Prevent retinal tissue damage
- Ameliorate visual recovery

## Team


### Research partner



### Industrial partner



### Medical partner



Project supported by Swiss Commission for Technology and Innovation CTI

### References:

- [1] S. Firański, J. Rivera, M. Clogenson, C. Baur, A. Bertholds, P. Llosas, S. Henein: Flexure-based multi-degrees-of freedom in-vivo force sensors for medical instruments. Euspen's 16th Conference & Exhibition, Nottingham, UK, May 2016
- [2] Optical force sensing element and microsurgical instruments, EP 12171195.6, Bertholds A., Llosas P. and Henein S., Patent Pending. Holder: Sensoptic SA (2012)
- [3] Optical measuring element having an integral structure, US12/919.621 EP 02255170, Bertholds A. & Llosas P., Henein S., Patent Holder: Sensoptic SA (Losone, CH) (2009).